



## Neo product series

### **Elegant and compact design**

The ergonomic and elegant design of the Neo dental treatment unit makes your practice more functional and convenient than ever before. The Neo's compact and small base provides for more working space for both the dentist and the assistant.

### **Hygienic, easy-to-clean surfaces**

Flush surfaces and seamless upholsteries make cleaning and disinfection easier and ensure good hygiene. The instrument bridge and suction head can be equipped with removable hygiene covers that are easy to disinfect.

### **Extensive movement range**

The extensive movement range of the Neo's arms ensures easy, comfortable operation and good work ergonomics. There are plenty of adjustment options to further ensure comfort. The large height adjustment range (45 to 95 cm) is designed to ensure a comfortable working position while standing or sitting. The chair can be rotated 90°, the suction head is adjustable from the 3 o'clock to 9 o'clock positions and suitable for both left- and right-handed use. Also, the cuspidor turns 90° sideways, providing more work space for the assistant.

### **Versatile and individualized adjustments**

Users can store individual settings for chair position and instrument speed ranges, among other settings. All user interfaces are clear and user-friendly. The Neo dental treatment system can be controlled with a wireless foot pedal, joystick or using keys on the user panels.

### **Everything within reach**

Instruments, trays, suction head, and cuspidor are all easy to reach and adjustable according to the needs of the dentist, the patient and the requirements of the current treatment. Thanks to the unit's thin headrest and backrest, the equipment never obstructs work and allows the dentist access to the treated area in an ergonomic and efficient position.